

CellMite® LVDT AC Excitation Dual Channel Intelligent Digital Signal Conditioner


Features:

- Dual Channel (M4338-200)
- 2-pt, mV/V Calibration
- 6-pt Linearization Calibration
- Scalable $\pm 10V$ Analog Output
- 16-Bit Analog Output
- Direct to PC with RS232 Data
- 24-Bit Resolution
- Multiple Filter Selections
- Storage for Three Calibrated Transducers
- Selectable Sensor Gains
- Remote Sense Excitation
- Addressable for Serial Network
- Tare, Peak, Valley.
- Unit Support: In, Cm, mm, %, Lb, Kg, N, mV/V, User Defined.
- TEDS-Tag™ Auto Identification
- Power Supplies Internally Isolated
- Graphical User Interfaces (GUI) and Software Drivers sold separately. Select from:
 - CellView Lite GUI (Cat. No 514569)
 - CellView Multi-Display GUI (Cat No. 514424)
 - Software Drivers (Cat No. 514864)

The Model 4338-200 Dual Channel CellMite® LVDT digital signal conditioner from Electro Standards Laboratories connects directly to two standard AC LVDT (Linear Variable Differential Transformer) transducers and simultaneously generates a serial output for connection to a PC and a $\pm 10V$ analog output. The CellMite® LVDT can be operated with a computer or simply setup by a computer and then used as a stand-alone intelligent LVDT digital signal conditioner with analog output.

The Model 4338-200 CellMite® LVDT has on board storage for up to three calibrated LVDT's, each with integrated 6-point calibration that can be used to linearize the LVDT outputs. Its multidrop RS232/RS485 serial port and simple command set allow for connection of multiple units in a serial network configuration. This compact unit is ideally suited for in-situ transducer conditioning and distributed process measurements.


Model 4338-200 features include nonvolatile memory for parameter and calibration storage, the ability to select between 3 stored LVDT calibrations, and automatic LVDT identification using TEDS-Tag™ technology. Solid state relay switches, multi-point and mV/V calibration, and the ability to locate LVDT core null position are incorporated into the unit. Internal generation of precision sine wave AC excitation voltages, remote sense excitation, and 24-bit internal resolution with 16-bit analog output are also standard.





Electro Standards Laboratories

ADVANCED SYSTEMS DESIGN & SERVICES



CellMite® 4338-200 LVDT Specifications:

Excitation:

Voltage: 3kHz, 3Vrms (std)
Sensors: LVDT's AC Bridges

Operation:

Internal Resolution: 24-bit.
Input Range: software adjustable from <10mV/V to >100mV/V
Conversion Rate: 60 per second single channel only
7.5/sec dual channel
Error: 0.05% +/-1 count
Precision Shunt: 60KΩ for bridge applications

Pushbuttons: Tare, Shunt

Indicator: Status LED

Outputs:

Analog: 16-bit, Scalable, ±10V.
Serial Data: Multi-drop RS232, RS485
Switch: Dual Solid State Relays

Power

Included Adapter: 5 VDC, 0.6 A (std), 110/240VAC, 60/50 Hz.
Optional Adapter: 12, 24, 48 VDC


Mechanical:

Size: 6" x 2.75" x 1.2"
Weight: 8.2 oz. (235 g)

CellView Lite GUI Software (Cat No. 514569)

Features:

- Turnkey data acquisition system
- Save data, calibration, and sensor test information to Excel spreadsheets
- Guides user through sensor calibration
- Guides user through adding/removing networked CellMite® and CellMite® LVDT units.
- Ability to select any networked M4338 unit.
- Ability to individually name M4338 units.
- Tare and reset peak and valley for the sensor.
- Save/Load the M4338 unit and sensor setup information.
- Set output data with user defined 5th order polynomial function.
- Setup a test to start and/or stop automatically
- Control M4338 solid state output relay switches.
- User programmable analog output voltage range.
- Units of In, Cm, mm, %, Lb, Kg, User Defined.



CellView Multi-Unit GUI Software (Cat No. 514424)

Features:

- Turnkey simultaneous data acquisition system for networked CellMite® and CellMite® LVDT units.
- Supports mix and match of CellMite® and CellMite® LVDT units on network.
- Ability to define 8 simultaneous data displays from multiple channels on multiple units.
- Save data, calibration, and sensor test information to Excel spreadsheets for all networked units.
- Guide user through sensor calibrations.
- Guides user through adding/removing networked CellMite® and CellMite® LVDT units.
- Ability to select and name any networked units.
- Tare and reset peak and valley for the sensors.
- Save/Load units and sensor setup information.
- Set output data with user defined 5th order polynomial function.
- Setup a test to start and/or stop automatically with user selectable trigger events and delayed triggering.
- Multi sample rate data acquisition via user definable logging profile.
- User programmable analog output voltage range for attached CellMite® and CellMite® LVDT units.
- Units of In, Cm, mm, %, Lb, Kg, mV/V, and User Defined.



36 Western Industrial Drive, Cranston, RI 02921
Tel: 401-943-1164 Fax: 401-946-5790

www.ElectroStandards.com
E-mail: eslab@ElectroStandards.com